



CLIMMATE SET-AW

Handles Even Sample Chunks

WATER ACTIVITY MEASUREMENT

QUICK FACTS

CLIMMATE SET-AW

Handheld Water Activity Meter

Handheld unit with big sample box to measure water activity of bigger amounts of sample or sample chunks. Can be used for raw material inspection or in production. Suitable to get a sneak peek what the water activity of the sample is about.

ADVANTAGES & BENEFIT

FOR BIGGER SAMPLE CHUNKS

To measure a sample in a representative way, bigger amount of sample is required sometimes. In this case, Novasina's ClimMate Set aw offers a comfortable solution. It allows the user to measure bigger quantity of sample or chunks of sample or even bulky stuff. This water activity meter does not feature the Novasina equilibrium detection, thus end of measurement has to be determined by the user himself and there is as well any implemented temperature-control. This meter is battery powered, so it can easily be used as portable version.

FEATURES

EASE OF USE	● ● ● ● ●
SIZE OF DISPLAY	● ● ● ● ●
DATA STORAGE	● ● ● ● ●
MEASUREMENT SPEED	● ● ● ● ●
MEASUREMENT ACCURACY	● ● ● ● ●
QUICK MODE	● ● ● ● ●
APPLICATION RANGE	● ● ● ● ●
CHEM. PROTECTION FILTER	● ● ● ● ●
THERMOSTATTED CHAMBER	● ● ● ● ●
PRICE	● ● ● ● ●
RE-USABLE STANDARDS	● ● ● ● ●
EASE OF MAINTENANCE	● ● ● ● ●
HUMIDITY EQUILIBRIUM DETECTION	● ● ● ● ●

novasina
The Art of Precision Measurement