

Disintegration Testers

6 Station (ED-6AL)

6
Station



Print




4 Station With Individual Timer
also available

Key Features

- 3 pair of baskets are controlled by independent timer and motor
- Snap click™ loading and swivel free basket movement
- Green illumination for better visibility of disintegration process
- Programmable temperature controller
- Magnetically coupled pump for water circulation in the bath
- Auto parking of the baskets at the top position at the end of the test
- Dual range timer
- Maintenance free smooth drive
- Basket assembly for Isotretinoin capsule analysis (optional)

Accessories



Disc handling
device (SAPO)



Magnetic
basket (SAPO)



Bolus basket with
3 tubes and guided discs



Bolus basket with
1 tube and guided disc



Basket cover