



Online Catalog

Pursuing a bright future with mixing technology

THINKY

Oh!

Key Factor for Achieving Maximum Performance of Planetary Centrifugal Mixers



Container / Adapter List for THINKY MIXER

“Are You Having Unsuccessful Results?”

Your containers or adapters could be the factor.

The proper **containers** and **adapters** are necessary for achieving 100% performance of a planetary centrifugal mixer.

THINKY provides original containers and adapters to fit the characteristics of the material.

We produce more than 150 custom-made adapters a year to meet customer needs.

Creating a whole new adapter is always challenging.

Our professional team considers the material characteristics, customer issues and the operating environment in order to design and supply you with custom-made adapters for your materials.

We are grateful for the frequent compliments from customers who appreciate the high quality of adapters made by THINKY.



Standard Adapters

● We provide various adapters for containers. Contact us for items not listed.



For cooling

250AD-100DSP-CO

Cooling adapter
for 100ml disposable cup



For cooling

500AD-300-CO

Cooling adapter
for CH300 container

Examples of custom made adapters

500AD-300CO-BHB

Cooling adapter for BHB 300ml anti-UV container

500AD-ASONE250CO

Cooling adapter for BOTTLES Wide Mouth 250ml (AS ONE)



100AD-100DSP

For 100 ml disposable cup



250AD-201

For 150 ml container



500AD-250

For 250 ml/300 ml container



250AD-100DSP-MP

For 100 ml disposable cup
(with lid)



250AD-200DSP-BLOW

For 200 ml disposable cup
(with lid)



250AD-100DSP

For 100 ml disposable cup



250AD-200DSP

For 200 ml disposable cup



500AD-500DSP3

For 500 ml disposable cup



250AD-205KAMI-CUP

For 205 ml paper cup



250AD-50SCR

For 50 ml vial glass bottle



250AD-200CH

For CH200 container



250AD-300CH

For CH300 container



500AD-ASONE500

For BOTTLES Wide Mouth
500ml (AS ONE)



250AD-20DT

For Dotite 20g



Heat-resistant Adapter

For 150 ml SUS container

Other standard adapters

250AD-20SCR

For 20 ml vial glass bottle

250AD-30SCR

For 30 ml vial glass bottle

Examples of custom-made adapters

250AD-GR10

For 10 ml vial PET bottle

250AD-GR30

For 30 ml vial PET bottle

250AD-BHB160

For BHB160 ml anti-UV
container

* All the items shown in the pictures above are sold separately.

Syringe Adapters

● We provide various adapters for syringes and barrels. Contact us for items not listed.



250AD-30S
For 20/30ml syringe



100AD-3S
For 3ml syringe



100AD-5S
For 5ml syringe



100AD-10S
For 10ml syringe



100AD-30S
For 20/30ml syringe



250AD-3S
For 3ml syringe



250AD-5S
For 5ml syringe



250AD-10S
For 10ml syringe



250AD-50S
For 50ml syringe



250AD-100BS
For 100ml barrel



500AD-170B
For 170ml barrel



250AD-M3S×6
For Musashi syringe
3ml□6

* The syringes/barrels shown in the pictures above are not included.

Other standard adapters

250AD-M5S×5
For Musashi syringe 5ml□5

250AD-M10S×4
For Musashi syringe 10ml□4

250AD-55S-EFD
For ex-EFD syringe 55ml

250AD-EFD3S×7
For ex-EFD syringe 3ml□7

250AD-EFD5S×5
For ex-EFD syringe 5ml□5

250AD-EFD10S×4
For ex-EFD syringe 10ml□4

Solder Paste Container Adapters

● We provide special container adapters for Solder Paste Mixer, SR-500.
(For a 500 g solder paste container that is produced in Japan, use 250AD-201 adapter on page 2.)

Examples of custom-made adapters

SP500AD-500LT
For LOCTITE/white-cap 500g

SP500AD-KESTER
For KESTER 500g


SP500AD-500LT-E
For LOCTITE/green-cap 500g

SP500AD-KESTER250
For KESTER 250g

SP500AD-HE
For HERAEUS 500g

Ointment Container Adapters

- We provide adapters for 12-125 ml straight-sided containers.
(For small volume mixing, the containers are convenient to use as disposable items.)
Contact us for items not listed.



250AD-NAN-U

For Umano UG 12 ml, 24 ml,
35 ml, 58 ml, 125 ml



50AD-NAN-U58

For Umano UG 58 ml



100AD-NAN-U

For Umano UG 12 ml, 24 ml,
35 ml

Example of Ointment Containers



Containers & Cups

● Disposable cups and plastic containers.



300 ml container



100 ml disposable cup



200 ml disposable cup



150 ml container



250 ml container



300 ml anti-UV container

* 160 ml anti-UV container is also available.



360 ml container [Type K] With keyways

* 360 ml container without keyways [Type F] is also available.



550 ml container



650 ml container



750 ml container



1000 ml container



4000 ml container



300 ml SUS
300 ml SUS container

Standard specifications

Without handle, without hole

Option

- (1) With handle, with hole
- (2) With handle, without hole
- (3) Without handle, with hole

* The picture above shows (2).

* Other types of containers (eg. container lid with a hole) are also available. Contact us for details.

Correspondence Table

○: Adapter required.

| container \ compatible model | AR-100 | ARE-310 ARE-250 AR-250 | ARE-400TWIN | ARE-500 AR-500 | ARV-310 ARV-310LED | ARV-5000 ARV-3000TWIN | SR-500 SP-500 | ARV-930TWIN |
|------------------------------|--------|------------------------------|-------------|-------------------|-----------------------|--------------------------|------------------|-------------|
| 100ml disposable cup | ○ | ○ | ○ | ○ | ○ | | | |
| 200ml disposable cup | | ○ | ○ | ○ | ○ | | | |
| 500ml disposable cup | | | | ○ | | | | ○ |
| 150ml container | | ○ | ○ | ○ | ○ | ○ | ○ | |
| 250ml container | | ● | ● | ○ | ● | ○ | ● | |
| 300ml container | | ● | ● | ○ | ● | ○ | ● | ○ |
| 300ml anti-UV container | | ● | ● | ○ | ● | ○ | ● | ○ |
| 360ml container | | | | ○ | | ○ | | ○ |
| 550ml container | | | | ● | | ○ | | ○ |
| 650ml container | | | | ● | | ○ | | |
| 750ml container | | | | | | ○ | | ● |
| 1000ml container | | | | | | ○ | | |
| 4000ml container | | | | | | ● | | |

Custom-made Adapters

THINKY designs more than 150 custom-made adapters every year. We create adapters especially suited for your containers. Contact us for more information.

Have you experienced the following?

Case 1

"We cannot find an adapter for our current container and have to transfer the material to a different container."



Custom-made adapters

for direct mixing and defoaming

Case 2

"Only a small amount needs to be processed but the container is too large for mixing small amounts."



Custom-made adapters

for mixing even 1ml of materials

Contact the address shown on the back cover of this brochure or our sales personnel.

Production flow of custom-made adapters



Contact

Contact the address shown on the back cover of this brochure or our sales personnel.



Interview and estimation

Sample containers are needed.



Meeting for detailed specifications



Order



Production



Completion

Example of custom-made adapters



The use of the most suitable containers and adapters results in optimum material processing. We consider the material characteristics, customer issues and the operating environment in order to design and supply you with custom-made adapters for your materials. Rapid response and attentive service are our strengths as a pure play company.

| Type | Model | Mixer | | | | | | Container | Remark |
|-------------------------|------------------------|--------|------------------------------|-----------------|-------------------|-----------------------|------------------------------------|--|----------------|
| | | AR-100 | ARE-310 ARE-250 AR-250 | ARE-400 TWIN | ARE-500 AR-500 | ARV-310 ARV-310LED | SR-500 SP-500 | | |
| For Standard Containers | 100AD-100DSP | ● | | | | | | 100 ml disposable cup | |
| | 250AD-100DSP | | ● | ● | ● | ● | ● | 100 ml disposable cup | |
| | 250AD-100DSP-MP | | ● | ● | ● | ● | ● | 100 ml disposable cup | |
| | 250AD-100DSP-CO | | ● | ● | ● | ● | ● | 100 ml disposable cup | cooling |
| | 250AD-200DSP | | ● | ● | ● | ● | ● | 200 ml disposable cup | |
| | 250AD-200DSP-BLOW | | ● | ● | ● | ● | ● | 200 ml disposable cup | |
| | 250AD-201 | | ● | ● | ● | ● | ● | 150 ml container 500 g Solder paste | |
| | 500AD-500DSP3 | | | | ● | | | 500 ml disposable cup | |
| | 500AD-250 | | | | ● | | | 250 ml container, 300 ml container | |
| | 250AD-20DT | | ● | ● | ● | ● | ● | DOTITE 20/30 g | |
| | 250AD-300SUS | | ● | ● | ● | ● | ● | 300 ml SUS | |
| | 250AD heat-resistant | | ● | ● | ● | ● | ● | 150 ml SUS | heat-resistant |
| | 250AD-205KAMI-CUP | | ● | ● | ● | ● | ● | 205 ml Paper cup | |
| | 250AD-GR10B | | ● | ● | ● | ● | ● | 10 ml PET bottle | |
| | 250AD-GR30B | | ● | ● | ● | ● | ● | 30 ml PET bottle | |
| | 250AD-20SCR | | ● | ● | ● | ● | ● | 20 ml glass bottle | |
| | 250AD-30SCR | | ● | ● | ● | ● | ● | 30 ml glass bottle | |
| | 250AD-50SCR | | ● | ● | ● | ● | ● | 50 ml glass bottle | |
| | 250AD-200CH | | ● | ● | ● | ● | ● | CH200 (TAKEMOTO) | |
| | 250AD-300CH | | ● | ● | ● | ● | ● | CH300 (TAKEMOTO) | |
| 250AD-BHB160 | | ● | ● | ● | ● | ● | BHB160 anti-UV | | |
| 500AD-ASONE500 | | | | | ● | | BOTTLES Wide Mouth 500 ml (AS ONE) | | |
| 500AD-300CO | | | | | ● | | CH300 (TAKEMOTO) | cooling | |
| 500AD-300CO-BHB | | | | | ● | | BHB300 anti-UV | cooling | |
| 500AD-ASONE250CO | | | | | ● | | BOTTLES Wide Mouth 250 ml (AS ONE) | cooling | |
| For Syringes & Barrels | 100AD-3S, 5S, 10S, 30S | ● | | | | | | 3ml, 5ml, 10ml, 20/30ml Syringe (ex-EFD, MUSASHI, IEI) | |
| | 250AD-3S, 5S, 10S, 30S | | ● | | ● | ● | ● | 3ml, 5ml, 10ml, 20/30ml Syringe (ex-EFD, MUSASHI, IEI) | |
| | 250AD-50S | | ● | ● | ● | (*) | (*) | 50 ml Syringe (MUSASHI, IEI) | |
| | 250AD-M3Sx6 | | ● | ● | ● | ● | ● | Musashi 3 ml Syringe x 6 (for NC-3E Syringe tip cap) | |
| | 250AD-M5Sx5 | | ● | ● | ● | ● | ● | Musashi 5 ml Syringe x 5 (for NC-3E Syringe tip cap) | |
| | 250AD-M10Sx4 | | ● | ● | ● | (*) | ● | Musashi 10 ml Syringe x 4 (for NC-3E Syringe tip cap) | |
| | 250AD-55S-EFD | | (*) | | | (*) | (*) | ex-EFD 55 ml Syringe | |
| | 250AD-EFD3Sx7 | | ● | ● | ● | ● | ● | ex-EFD 3 ml Syringe x 7 | |
| | 250AD-EFD5Sx5 | | ● | ● | ● | ● | ● | ex-EFD 5 ml Syringe x 5 | |
| | 250AD-EFD10Sx4 | | ● | ● | ● | (*) | ● | ex-EFD 10 ml Syringe x 4 | |
| | 250AD-100BS | | ● | ● | ● | ● | ● | 100 ml Barrel (ex-EFD, MUSASHI, PPG) | |
| | 500AD-170B | | | | | ● | | 170 ml Barrel (ex-EFD, MUSASHI, PPG) | |
| For Solder Paste | SP500AD-500LT | | | | | | ● | * LT (Wc) 500 g | |
| | SP500AD-500LT-E | | | | | | ● | * LT (Gc) , Cs, SH, SM, MIC 500 g | |
| | SP500AD-HE | | | | | | ● | * HE 500g | |
| | SP500AD-KESTER | | | | | | ● | * KS 500g | |
| | SP500AD-KESTER250 | | | | | | ● | * KS 250g | |
| For Ointment Containers | 50AD-NAN-U58 | ● | | | | | | UG 58ml | |
| | 100AD-NAN-U | ● | | | | | | UG 12ml, 24ml, 35ml | |
| | 250AD-NAN-U | | ● | ● | ● | ● | | UG 12ml, 24ml, 35ml, 58ml, 125ml | |

* Solder Paste Maker LT (Gc): LOCTITE (Green cap) LT (Wc): LOCTITE (White cap) HE: HERAEUS KS: KESTER Cs: Cookson SH: SHENMAO SM: SM泰克 MIC: MICROSDLER
 ■ Product appearance/specifications may change without notice.

For requests concerning demonstrations and evaluation testing, please contact
THINKY CORPORATION

or the sales agent below

For the latest information about products and exhibitions, visit:

<https://www.thinkymixer.com/en-gl/>

THINKY CORPORATION

Headquarters: 2-16-2 Sotokanda, Chiyoda-ku, Tokyo 101-0021
 Phone : +81-3-5207-2713 Fax : +81-3-5289-3281

THINKY USA Inc. : 23151 Verdugo Drive, Suite 112 Laguna Hills, CA92653, USA

THINKY CHINA: East building, HaiAn Kafunuo Mansion, Shennan road,
 Qianhai road, Nanshan district, Shenzhen

Sales agent